

KEFĂ

 **phenomenex**[®]
*...breaking with tradition*SM

10% POPUST

Ponuda važi do 30.09.2024.



roQTM QuEChERS Kits

An Easier
QuEChERS Solution
for Diverse Matrices

 **phenomenex**[®]
*...breaking with tradition*SM

Learn more at
www.phenomenex.com/roQ




What is roQ™ QuEChERS?

QuEChERS is a descriptive acronym that stands for Quick, Easy, Cheap, Effective, Rugged, and Safe. Improved with you in mind, the unique design of the roQ QuEChERS kits eliminates common problems seen with current QuEChERS kits. By incorporating features such as trays to hold centrifuge tubes, easy pour salt packets, and stand alone centrifuge tubes, roQ QuEChERS Kits make the QuEChERS procedure even easier!

Stay Organized with Built-In Test Tube Racks

roQ boxes are designed with removable racks that can be used as test tube racks




Easily Weigh Out Samples with Stand Alone Centrifuge Tubes

Simultaneously add sample while weighing, no transfer steps required



Cleaner Extracts with Low Extractable Tubes

Reduce background noise in chromatograms



15 mL roQ centrifuge tube extracted with 1 % Acetic acid in Acetonitrile and then with Toluene.


No Salt Spills with Easy Pour Salt Packets

Our easy open, easy pour salt packets have a narrow spout which makes salt addition easy



Avoid Leaky Tubes with No Leak Caps

We guarantee our tubes won't leak!




How does it work?

The QuEChERS technique addresses shortcomings of traditional sample preparation methods, such as long extraction procedures and the use of hazardous solvents, and radically simplifies multi-residue pesticide analysis in food or environmental samples.

The QuEChERS technique consists of two steps; Extraction and Dispersive SPE (dSPE)


Step 1: Extraction

Pesticides and analytes of interest must first be extracted from the sample. This process relies on the combination of organic solvent and various salts to partition the analytes from samples into an organic layer (typically acetonitrile).




Step 2: Dispersive SPE (dSPE)

An aliquot of the organic layer from the extraction step is subjected to further clean up by dispersive SPE. This step selectively removes unwanted interferences such as lipids and pigments.




Add
supernatant from extraction procedure into a roQ dSPE tube.

Avoid Leaky Tubes
Our tubes are designed to seal completely




Shake
dSPE tube for 30 seconds.



Centrifuge
dSPE tube for 5 minutes.

Cleaner Extracts
roQ kits contain low extractable centrifuge tubes reducing background noise in chromatograms



Analyze
supernatant by GC or HPLC.

Step 1

roQ™ Extraction Kits

QuEChERS can be performed by following 3 different methods: the AOAC 2007.01 Method, the EN 15662 Method, or the Original Non-Buffered Method.

Use the chart below to select the most appropriate kits for your work and sample size/type.

AOAC 2007.01 Method

Salts used:

Magnesium Sulfate (MgSO₄)

Induces phase separation between water content in sample and acetonitrile layer

Sodium Acetate (NaOAc)

Buffers the sample to stabilize pH

Original Non-Buffered Method

Salts used:

Magnesium Sulfate (MgSO₄)

Induces phase separation between water content in sample and acetonitrile layer

Sodium Chloride (NaCl)

Induces phase separation between water content in sample and acetonitrile layer

EN 15662 Method

Salts used:

Magnesium Sulfate (MgSO₄)

Induces phase separation between water content in sample and acetonitrile layer

Sodium Chloride (NaCl)

Induces phase separation between water content in sample and acetonitrile layer

Sodium Citrate Tribasic Dihydrate (SCTD)

Buffers the sample to stabilize pH

Sodium Citrate Dibasic Sesquihydrate (SCDS)

Buffers the sample to stabilize pH

Sample Size

15 g

Part No.: **KS0-8911**

15 g

Part No.: **KS0-8912**

to be used with AOAC 2007.01 dSPE procedure

or

10 g

Part No.: **KS0-8910**

to be used with EN 15662 dSPE procedure

10 g

Part No.: **KS0-8909**

All roQ Extraction Kits include 50 easy-pour salt packets and fifty 50 mL stand-alone centrifuge tubes.

Step 2

Select Your roQ™ dSPE Kit

Fruit and Vegetable Type

Method	Sample Size	General		Fats and Waxes		Pigmented		Highly Pigmented		Pigments and Fats	
		Contents	Part No.	Contents	Part No.	Contents	Part No.	Contents	Part No.	Contents	Part No.
AOAC 2007.01	1 mL	150 mg MgSO ₄ 50 mg PSA	KS0-9511	150 mg MgSO ₄ 50 mg PSA 50 mg C18E	KS0-9512	150 mg MgSO ₄ 50 mg PSA 50 mg GCB	KS0-9513	—	—	150 mg MgSO ₄ 50 mg PSA 50 mg GCB 50 mg C18E	KS0-9514
	8 mL	1200 mg MgSO ₄ 400 mg PSA	KS0-9515	1200 mg MgSO ₄ 400 mg PSA 400 mg C18E	KS0-9516	1200 mg MgSO ₄ 400 mg PSA 400 mg GCB	KS0-9517	—	—	1200 mg MgSO ₄ 400 mg PSA 400 mg GCB 400 mg C18E	KS0-9518
EN 15662	1 mL	150 mg MgSO ₄ 25 mg PSA	KS0-9503	150 mg MgSO ₄ 25 mg PSA 25 mg C18E	KS0-9504	150 mg MgSO ₄ 25 mg PSA 2.5 mg GCB	KS0-9505	150 mg MgSO ₄ 25 mg PSA 7.5 mg GCB	KS0-9506	—	—
	6 mL	900 mg MgSO ₄ 150 mg PSA	KS0-9507	900 mg MgSO ₄ 150 mg PSA 150 mg C18E	KS0-9508	900 mg MgSO ₄ 150 mg PSA 15 mg GCB	KS0-9509	900 mg MgSO ₄ 150 mg PSA 45 mg GCB	KS0-9510	—	—

Salts and sorbents used in roQ dSPE kits

Magnesium Sulfate (MgSO₄)	Removes excess water from sample
Primary/Secondary Amine (PSA)	Removes organic acids, fatty acids, sugars, and anthocyanin pigments from sample
Endcapped C18 Sorbent (C18E)	Removes fats, sterols, and other non-polar interferences from sample
Graphitized Carbon Black (GCB)	Removes pigments from sample NOT FOR USE WITH PLANAR PESTICIDES

All roQ dSPE Kits include pre-weighed sorbents/salts inside 2 mL or 15 mL centrifuge tubes.

Extraction of Pesticide Residues from Kale and Grapes Using roQ™ EN 15662 QuEChERS Kit

LC-MS Analysis

Over 200 pesticides were screened at concentration ranges between 0.01 ppm to 1 ppm with the majority of analytes having a recovery range of 70-130%. The method produced excellent selectivity and reproducibility for the earlier eluting polar pesticides owing to the use of QuEChERS for removing sample matrix interferences and the biphenyl bonded-phase chemistry of the Kinetex® Biphenyl LC column for the chromatographic separation across the entire range of pesticides.

QuEChERS Procedure (EN 15662 Method)

Kale and grapes samples were prepared at four concentrations: no spike (0.0 ppm), low spike (0.1 ppm), medium spike (0.5 ppm), and high spike (1 ppm)

Pretreatment: Kale and grapes were frozen at ~-80°C overnight

Step 1: Extraction

1. Homogenize sample using a blender or similar apparatus.
2. Weigh 10 g of homogenized sample into a clean 50 mL tube (provided in roQ Extraction Kits).
3. Add 10 mL of Acetonitrile containing internal standard.
4. Dispense contents of the included QuEChERS salt packet into the 50 mL tube containing homogenized sample.
5. Shake vigorously by hand for 1 minute.
6. Centrifuge for 5 minutes @ 4000 rpm, making sure that the solid material is at the bottom of the tube and a liquid layer forms on top of the solid material.

Step 2: Dispersive SPE (dSPE)


1. Transfer 6 mL of supernatant from Step 6 of the extraction process into the 15 mL tube containing the QuEChERS dSPE sorbents.
2. Shake vigorously by hand for 30 seconds.
3. Centrifuge for 5 minutes @ 4000 rpm to separate solid material from the liquid layer.
4. Transfer desired supernatant to an autosampler vial.
5. Dilute samples 1:10 in mobile phase A prior to injection or LC-MS analysis.

To read full technical note, go to
www.phenomenex.com and


Search:



Standard Chromatogram (0.5 ppm standard)




Kale Sample (1 ppm spike)




LC-MS/MS Conditions

Column: Kinetex® 5 µm Biphenyl	Flow Rate: 0.5 mL/min
Dimensions: 50 x 4.6 mm	Inj. Volume: 20 µL
Part No.: 00B-4627-E0	Temperature: 35 °C
SecurityGuard™ ULTRA: AJ0-9207	Detection: MS/MS (ESI+), scheduled MRM
Mobile Phase: A: 10mM Ammonium Formate in Water	Instrument: SCIEX 4000 QTRAP®
B: Methanol	
Gradient:	
Time (min) % B	
0 2	
1 2	
10 100	
13 100	
13.1 2	
16 2	

Sample: For full list of 210 pesticides, go to www.phenomenex.com and Search: 



Grapes Sample (0.5 ppm spike)



Mycotoxins from Corn Meal Products Using roQ™ EN 15662 QuEChERS Kit LC-MS/MS

The use of roQ QuEChERS and Kinetex® XB-C18 Core-Shell Technology LC columns deliver a rapid and simple approach for mycotoxin screening from corn products.

QuEChERS Procedure (EN 15662 Method)

Step 1: Extraction from Ground Corn

1. Homogenize sample using a blender or similar apparatus
2. Weigh and transfer 5g of ground corn meal to a 50mL roQ QuEChERS extraction tube
3. Add 10mL of water and 10mL of acetonitrile with 1.0% formic acid
4. Dispense contents of the included roQ QuEChERS extraction packet (KS0-8909) into the 50mL tube containing homogenized sample
5. Shake vigorously by hand for 1 minute
6. Centrifuge for 5 minutes @ 4000 rpm, making sure that the solid material is at the bottom of the tube and a liquid layer forms on top of the solid material

Step 2: Clean up using dispersive Solid Phase Extraction (dSPE)

1. Transfer the supernatant from Step 6 of the extraction process into a roQ QuEChERS 15mL centrifuge tube containing 900mg MgSO₄ and 150mg PSA (KS0-9507)
2. Shake vigorously by hand for 30 seconds
3. Centrifuge for 5 minutes at 4000 rpm to separate solid material from the liquid layer
4. Transfer the supernatant to a vessel for evaporation


LC-MS/MS Conditions

Column: Kinetex® 2.6µm XB-C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4496-AN
SecurityGuard Cartridge: AJO-8782
Mobile Phase: A: 5 mM Ammonium acetate with 0.5% Acetic acid
 B: 5 mM Ammonium acetate in Methanol with 0.5% Acetic acid


Gradient: Time (min)	% B
0	5
2	5
5	80
5.2	98
8	98

Flow Rate: 0.45 mL/min
Temperature: Ambient
Injection Volume: 25 µL
Detection: Tandem Mass Spec (MS/MS)
Instrument: SCIEX API 5000™

Positive Mode



Negative Mode



Free Method Development Support and Guidance



Did you know?

You have a Phenomenex technical support team dedicated to supporting your lab's needs! We collaborate with chemists every day to improve productivity and results while working within the bounds of your established methods.

Email: info@phenomenex.com



Determination of PFASs in Sediments Using roQ™ AOAC 2007.01 QuEChERS Kits

LC-MS Analysis

In order to extract Perfluoroalkyl substances (PFASs) from marine and freshwater sediment, QuEChERS was introduced to save time and offer reliable results.

QuEChERS Procedure (AOAC 2007.01 Method)

Step 1: Extraction from Sediments

1. Weigh 2.0 g of suitably dried sediment into a polypropylene container and spike with isotopically-labeled internal standards. PPCPs, Steroids, and Pyrethroids can be extracted concurrently with this method by adding the appropriate internal standard and spiking solutions to the samples and QCs.
2. Add 10 mL deionized water and vortex. Add 10 mL acidified acetonitrile (1 % Acetic acid) to the slurry and vortex.
3. Add the extraction salts (1.5 g Sodium acetate and 2 g MgSO₄) to the sample and vortex for 1 minute.
4. Centrifuge the samples for 5 minutes at 4000 rpm.
5. Place the samples in a rack and freeze at -20 °C for 30-60 minutes. This freezing step allows for easier extraction of the supernatant.

Step 2: Clean up Using dispersive Solid Phase Extraction (dSPE)


1. Transfer 8-9 mL of the acetonitrile supernatant into a roQ QuEChERS PSA/C18 dSPE clean-up tube (Part no. KS0-9516)
2. Vortex for one minute.
3. Centrifuge the dSPE tubes for 10 minutes at 3000 rpm.
4. Place an aliquot of the extract in a polypropylene HPLC vial and dilute 1:1 with deionized water. The sample is now ready for analysis.

LC-MS/MS Conditions

Column:	Gemini® 3 µm C18	
Dimensions:	100 x 3.0 mm	
Part No.:	00D-4439-Y0	
Inline Filter:	Phenomenex KrudKatcher™ Ultra (AFO-8497)	
Delay Column:	Luna® 5 µm C18(2) 30 x 2.0 mm	
Part No.:	00A-4252-B0	
Mobile Phase:	A: 20 mM Ammonium acetate in Water B: Methanol	
Gradient:	Time (min)	% B
	0	10
	1.5	65
	8	95
	8.1	99
	12	99
Flow Rate:	0.6 mL/min	
Injection:	90 µL	
Temperature:	40 °C	
Detector:	SCIEX 5500 QTRAP®	
Detection:	MS/MS ESI Negative (sMRM)	
Samples:	1. PFBA	11. PFNA
	2. PFPeA	12. PFOSA
	3. PFBS	13. PFNS
	4. PFHxA	14. PFDA
	5. PFPS	15. PFDS
	6. PFHxS	16. PFUdA
	7. PFHpA	17. PFDoA
	8. PFHpS	18. PFTrDA
	9. PFOA	19. PFTeDA
	10. PFOS	



Sediment Spiked with 1 ng/g



Method Performance Data for Sediments Spiked at 1 ng/g of the Target Analytes (n=4)

Compound	Average % Recovery	% RSD
PFBA	91.7	0.76
PFPeA	86.3	6
PFHxA	89.4	1.2
PFHpA	93.1	2.9
PFOA	98.3	1.5
PFNA	93	1.6
PFDA	87.7	4.5
PFUdA	92.3	2.1
PFDoA	92.5	4.1
PFTrDA	88.2	2.1
PFTeDA	87.6	2.1
PFBS	86.3	2.1
PFPeS	96.2	3.2
PFHxS	81.3	5
PFHpS	92.3	2.6
PFOS	92.1	2.6
PFOSA	104.5	6.3
PFNS	89.8	6.8
PFDS	87.3	6.7

Learn more about roQ QuEChERS

Visit: www.phenomenex.com/roQ



Acknowledgement

Special thanks to Syljohn Estil and the Sanitation Districts of Los Angeles County - San Jose Water Quality Laboratory for contributing this method.

Analysis Tools You Can Depend On

Achieve resolution of multi-residues with Phenomenex HPLC, UHPLC, and GC columns.

UHPLC/HPLC

Luna Omega High Performance Silica-based HPLC Columns



Enhanced with 20 years of technology, innovation, and experience, Luna Omega columns build upon the Luna legacy to now provide incredible UHPLC/HPLC performance and selectivity. With Luna Omega columns you gain:

- Industry benchmark for LC columns
- Wide pH stability for long column lifetime
- Extensive method validation document for proven reproducibility



Visit www.phenomenex.com/luna for more information.

UHPLC/ HPLC

Kinetex Core-Shell LC Columns



Leverage the power of Kinetex to achieve ultra-high performance on any LC system.

- Increase resolution, throughput, and sensitivity
- Decrease solvent consumption
- Save time and money
- Transfer methods anywhere



Visit www.phenomenex.com/kinetex for more information.

GC

Zebbron High Performance ZB-MultiResidue™ Columns



ZB-MultiResidue-1 & -2 columns are optimized for pesticides, herbicides, and insecticide analysis. These columns have an orthogonal selectivity that allows them to be used together in dual column confirmation analysis.

- Low activity, decreased breakdown of sensitive pesticides such as DDT
- MS Certified phases provide low bleed performance for pesticide confirmation by MS
- Long column lifetimes



Visit www.phenomenex.com/zebron for more information.



Ordering Info

roQ Extraction Kits

Extraction kits contain fifty easy-pour salt packets and fifty 50 mL stand-alone centrifuge tubes

Description	Unit	Part No.
AOAC 2007.01 Method Extraction Kits		
6.0 g MgSO ₄ , 1.5 g NaOAc	50/pk	KSO-8911
EN 15662 Method Extraction Kits		
4.0 g MgSO ₄ , 1.0 g NaCl, 1.0 g SCTD, 0.5 g SCDS	50/pk	KSO-8909
Original Non-buffered Method Extraction Kits		
4.0 g MgSO ₄ , 1.0 g NaCl	50/pk	KSO-8910
6.0 g MgSO ₄ , 1.5 g NaCl	50/pk	KSO-8912

roQ dSPE Kits

dSPE kits contain pre-weighed sorbents/salts inside 2 mL or 15 mL centrifuge tubes

Description	Unit	Part No.
2 mL dSPE Kits		
150 mg MgSO ₄ , 25 mg PSA, 25 mg C18E	100/pk	KSO-9504
150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB	100/pk	KSO-9505
150 mg MgSO ₄ , 25 mg PSA, 7.5 mg GCB	100/pk	KSO-9506
150 mg MgSO ₄ , 25 mg PSA	100/pk	KSO-9503
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18E, 50 mg GCB	100/pk	KSO-9514
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18E	100/pk	KSO-9512
150 mg MgSO ₄ , 50 mg PSA, 50 mg GCB	100/pk	KSO-9513
150 mg MgSO ₄ , 50 mg PSA	100/pk	KSO-9511
15 mL dSPE Kits		
900 mg MgSO ₄ , 150 mg PSA, 150 mg C18E	100/pk	KSO-9508
900 mg MgSO ₄ , 150 mg PSA, 15 mg GCB	100/pk	KSO-9509
900 mg MgSO ₄ , 150 mg PSA, 45 mg GCB	100/pk	KSO-9510
900 mg MgSO ₄ , 150 mg PSA	100/pk	KSO-9507
1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18E, 400 mg GCB	100/pk	KSO-9518
1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18E	100/pk	KSO-9516
1200 mg MgSO ₄ , 400 mg PSA, 400 mg GCB	100/pk	KSO-9517
1200 mg MgSO ₄ , 400 mg PSA	100/pk	KSO-9515

Bulk roQ QuEChERS Sorbents

Phase	10g	100g
C18-E	—	04G-4348
GCB (Graphitized Carbon Black)	04D-4615	04G-4615
PSA	—	04G-4610

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna, Gemini, and Kinetex are registered trademarks of Phenomenex. MultiResidue, KrudKatcher, SecurityGuard, Be-Happy, roQ, and Zebtron are trademarks of Phenomenex. API 5000 is a trademark and QTRAP is a registered trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

Disclaimer

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2018 Phenomenex, Inc. All rights reserved.

roQ Extraction Salt Packets

Salt packets only. Centrifuge tubes not included.

Description	Unit	Part No.
AOAC 2007.01 Method Extraction Packets		
6.0 g MgSO ₄ , 1.5 g NaOAc	50/pk	KSO-9043
EN 15662 Method Extraction Packets		
4.0 g MgSO ₄ , 1.0 g NaCl, 1.0 g SCTD, 0.5 g SCDS	50/pk	KSO-9041
Original Non-Buffered Method Extraction Packets		
4.0 g MgSO ₄ , 1.0 g NaCl	50/pk	KSO-9042
6.0 g MgSO ₄ , 1.5 g NaCl	50/pk	KSO-9044

We're here to help!

Contact your Sample Preparation Specialist
By email: Support@Phenomenex.com



BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

**All other countries
Corporate Office USA** 

t: +1 (310) 212-0555
info@phenomenex.com



www.phenomenex.com


Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

MYCOTOXINS

Analysis

Using QuEChERS, SPE & LC/MS/MS

- *Selective*
- *Rapid*
- *Repeatable*



 **phenomenex**[®]
...breaking with traditionSM



www.phenomenex.com/Food

Rapid and Selective Multi-Toxin Screening

Column choices for enhanced retention, selectivity and speed

Column selectivity choices such as Kinetex XB-C18 and Kinetex Biphenyl Core-Shell Technology columns can provide you with a wide range of retention capabilities for Mycotoxin screening.

LC/MS/MS Conditions

Column: Kinetex® 2.6 µm XB-C18 100 Å

Dimensions: 50 x 2.1 mm

Part No.: 00B-4496-AN

Mobile Phase: A: Water with 5 mM Ammonium acetate and 0.5% Acetic acid

B: Methanol with 5 mM Ammonium acetate and 0.5% Acetic acid

Gradient Time (min)	% B
0	2
2	2
5	80
5.2	98
8	98

Flow Rate: 450 µL/min


Temperature: Ambient

Detection: Tandem Mass Spectrometer (MS/MS) (550 °C)


Detector: AB SCIEX API 5500™

- Sample:**
- 15-Acetyldeoxynivalenol (15-AcDON)
 - 3-Acetyldeoxynivalenol (3-AcDON)
 - Aflatoxin B1 (AFB1)
 - Aflatoxin B2 (AFB2)
 - Aflatoxin G1 (AFG1)
 - Aflatoxin G2 (AFG2)
 - Deoxynivalenol (DON)
 - Diacetoxyscirpenol (DAS)
 - Fumonisin B1 (FB1)
 - Fumonisin B2 (FB2)
 - Fusarenone-X (FUS X)
 - HT-2 Toxin
 - Monoacetoxyscirpenol (MAS)
 - Nivalenol (NIV)
 - Ochratoxin A (OTA)
 - T-2 Toxin
 - Zearalenone (ZEA)

Negative Polarity



Positive Polarity



TIC of all analytes with negative and positive fast polarity switching.

App ID 20874



LC/MS/MS Conditions


Column: Kinetex 2.6µm Biphenyl
Dimensions: 50 x 2.1 mm
Part No.: 00B-4622-AN
Mobile Phase: A: 5 mM Ammonium acetate with 0.1% Acetic acid
 B: Methanol with 5 mM Ammonium acetate with 0.1% Acetic acid

Gradient:	Time (min)	% B
	0	2
	2	2
	5	80
	5.2	98
	8	98
	8.01	2
	11	2

Flow Rate: 0.45 mL/min
Temperature: Ambient
Detection: MS/MS (AB SCIEX API 4000™)

Sample:

- 15-Acetyldeoxynivalenol
- DAS
- FB1
- HT2 Toxin
- FB2
- T2 Toxin
- Aflatoxin M1
- Aflatoxin G2
- Ochratoxin A
- Aflatoxin G1
- Aflatoxin B2
- Aflatoxin B1



Expand your Selectivity Options with Kinetex Core-Shell Technology
www.phenomenex.com/Kinetex

Easy Extraction and Fast Screening from Corn Meal Products

The use of roQ QuEChERS and Kinetex XB-C18 Core-Shell Technology columns deliver a rapid and simple approach for Mycotoxin screening from corn products.

roQ
QuEChERS Kits



Sample Preparation

QuEChERS

Extraction from Ground Corn

1. Homogenize sample using a blender or similar apparatus
2. Weigh and transfer 5 g of ground corn-meal to a 50 mL roQ QuEChERS extraction tube
3. Add 10 mL of water and 10 mL of acetonitrile with 1.0 % formic acid
4. Dispense contents of the included roQ QuEChERS extraction packet (KS0-8909) into the 50 mL tube containing homogenized sample
5. Shake vigorously by hand for 1 minute
6. Centrifuge for 5 minutes @ 4000 rpm, making sure that the solid material is at the bottom of the tube and a liquid layer forms on top of the solid material

Quick and Easy
Procedure

Clean up using dispersive Solid Phase Extraction (dSPE)


1. Transfer the supernatant from Step 6 of the extraction process into a roQ QuEChERS 15 mL centrifuge tube containing 900 mg $MgSO_4$ and 150 mg PSA (KS0-8924)
2. Shake vigorously by hand for 30 seconds
3. Centrifuge for 5 minutes at 4000 rpm to separate solid material from the liquid layer
4. Transfer the supernatant to a vessel for evaporation

roQ QuEChERS Simplifies your Sample Preparation
visit www.phenomenex.com/roQ for details.

LC/MS/MS Conditions

Column:	Kinetex® 2.6 µm XB-C18												
Dimensions:	50 x 2.1 mm												
Part No.:	00B-4496-AN												
SecurityGuard Cartridge:	AJ0-8782												
Mobile Phase:	A: 5 mM Ammonium acetate with 0.5 % Acetic acid B: 5 mM Ammonium acetate in Methanol with 0.5 % Acetic acid												
Flow Rate:	0.45 mL/min												
Gradient:	<table><thead><tr><th>Time (min)</th><th>% B</th></tr></thead><tbody><tr><td>0.0</td><td>5</td></tr><tr><td>2.0</td><td>5</td></tr><tr><td>5.0</td><td>80</td></tr><tr><td>5.2</td><td>98</td></tr><tr><td>8.0</td><td>98</td></tr></tbody></table>	Time (min)	% B	0.0	5	2.0	5	5.0	80	5.2	98	8.0	98
Time (min)	% B												
0.0	5												
2.0	5												
5.0	80												
5.2	98												
8.0	98												
Injection Volume:	25 µL												
Temperature:	Ambient												
Detection:	API 5000 (AB SCIEX) Tandem Mass Spec (MS/MS)												
System:	Agilent® 1200SL LC system (Agilent Technologies, Palo Alto, CA, USA), equipped with a binary pump autosampler												

Figure 1. Sample Extract Ion Chromatogram
Positive Mode




- Sample:**
1. 15-Acetyldeoxynivalenol (15-AcDON)
 2. Aflatoxin G2 (AFG2)
 3. Aflatoxin G1 (AFG1)
 4. Unknown Interference
 5. Aflatoxin B2 (AFB2)
 6. Aflatoxin B1 (AFB1)
 7. Diacetoxyscirpenol (DAS)
 8. HT2 Toxin
 9. Fumonisin B1 (FB1)
 10. T-2 Toxin
 11. Ochratoxin A (OTA)
 12. Fumonisin B2 (FB2)

Table 1. MRM Transitions & Retention Times for Mycotoxins

Analyte	Q1	Q3	Retention Time (min)	Mode	% Recovery
Patulin	152.9	108.9/81.0	0.65	-ve	73
Nivalenol (NIV)	371	281.0/59.0	0.88	-ve	87
Deoxynivalenol (DON)	355	295.0/58.6	1.57	-ve	67
Fusarenone-X	413	59.1/353.1	3.76	-ve	93
3-Acetyldeoxynivalenol (3-AcDON)	397	59.0/307.0	4.34	-ve	96
Zearalenone (ZEA)	317	175.0/272.9	5.64	-ve	91
Verruculogen	510.1	166.2/305.9	5.96	-ve	90
Diacetoxyscirpenol (DAS)	384	307.0/349.0	5.05	+ve	79
15-Acetyldeoxynivalenol (15-AcDON)	339	321.0/137.0	4.35	+ve	99
Aflatoxin G2 (AFG2)	331	285.0/245.0	4.64	+ve	8
Aflatoxin G1 (AFG1)	329	243.0/311.0	4.73	+ve	14
Aflatoxin B2 (AFB2)	315.1	287.0/259.1	4.84	+ve	118
Aflatoxin B1 (AFB1)	313	285.1/257.1	4.94	+ve	91
HT2 Toxin	447.1	345.0/285.1	5.36	+ve	80
T2 Toxin	489	387.0/327.2	5.60	+ve	66
Fumonisin B1 (FB1)	722.2	352.2/528.2	5.81	+ve	96
Ochratoxin A (OTA)	404.1	239.0/386.1	5.60	+ve	83
Fumonisin B2 (FB2)	706.2	336.2/512.5	5.81	+ve	81

Figure 2. Sample Extract Ion Chromatogram
Negative Mode



- Sample:**
1. Patulin
 2. Nivalenol (NIV)
 3. Deoxynivalenol (DON)
 4. Fusarenone-X
 5. 3-Acetyldeoxynivalenol (3-AcDON)
 6. Zearalenone (ZEA)
 7. Verruculogen



Boost your Analysis Speeds with Kinetex Core-Shell Technology
See How www.phenomenex.com/kinetex

Aflatoxins, Ochratoxin A, Trichothecenes, Zearalenone in Cereal Based Foods

Strata-X solid phase extraction followed by LC/MS/MS offers reliable extraction and identification of key Mycotoxins in cereal based foods.




LC/MS/MS Conditions

Column: Kinetex® 2.6 µm C18 100Å
Dimensions: 100 x 2.1 mm
Part No.: 00D-4462-AN
Mobile Phase: A: Water with 0.5% Acetic acid and 1 mM Ammonium acetate
 B: Methanol with 0.5% Acetic acid and 1 mM Ammonium acetate
Gradient:

Time (min)	% B
0	10
4	40
31	60
36	60

Flow Rate: 200 µL/min
Temperature: 40 °C
Detection: Tandem Mass Spectrometer (MS/MS) (350 °C)

Analyte Name	RT (min)
Deoxynivalenol	4.8
Aflatoxin G2	12.8
Aflatoxin G1	13.6
Aflatoxin B2	14.6
Aflatoxin B1	15.8
HT-2 Toxin	20.1
T-2 Toxin	25.3
Zearalenone	29.0
Ochratoxin A	31.4



Total ion chromatogram (sum of MRM transitions) of a wheat based crisp bread sample extract spiked with:
 750 µg/kg DON; 1 µg/kg AFG2, AFB2; 3 µg/kg AFG1; 5 µg/kg AFB1; 50 µg/kg HT-2, T-2; 75 µg/kg ZEA; 3 µg/kg OTA

Author reference:
 Michele Suman
 michele.suman@barilla.com
 Barilla Food Research Labs,
 via Mantova 166, 43100 Parma, Italy





Sample Preparation

Solid Phase Extraction

Strata-X 60 mg/3 mL

Part No.: 8B-S100-UBJ
Condition: 2 mL Methanol
Equilibrate: 2 mL Methanol/Water (10:90)
Load: Pretreated sample
Wash: 1 mL Methanol/Water (20/80, v/v)
Elute: 1 mL of Methanol



Chromatograms (MRM transitions) of flour naturally contaminated with:
Top: 898 µg/kg DON; 9 µg/kg HT2; 2 µg/kg T2; 6 µg/kg OTA (soft wheat)
Bottom: 5.3 µg/kg DON; 5.4 µg/kg HT2; 1.5 µg/kg T2; 0.8 µg/kg OTA (durum wheat)

Veronica M. T. Lattanzio
 veronica.lattanzio@ispa.cnr.it
 Institute of Sciences of Food Production (ISPA),
 National Research Council of Italy (CNR),
 via Amendola 122/0, 70126 Bari, Italy



Ordering Information

Kinetex®

Get the most out of every HPLC/UHPLC analysis. Kinetex Core-Shell Technology produces increased efficiencies over traditional fully porous columns, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity.

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges [†]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
Biphenyl	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJO-9207
XB-C18	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJO-8768
C18	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJO-8768

for 4.6 mm ID



2.6 µm Minibore Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
Biphenyl	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJO-9209
XB-C18	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
C18	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
F5	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJO-9322

for 2.1 mm ID

[†]SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

Strata™-X

Versatile polymeric Solid Phase Extraction sorbent for effective retention of any analyte of interest.

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	30 mg	8B-S100-TAK**	1 mL (100/box)	
	30 mg	8B-S100-TBJ	3 mL (50/box)	
	60 mg	8B-S100-UBJ**	3 mL (50/box)	
	100 mg	8B-S100-EBJ	3 mL (50/box)	
	100 mg	8B-S100-ECH	6 mL (30/box)	
	200 mg	8B-S100-FBJ	3 mL (50/box)	
	200 mg	8B-S100-FCH	6 mL (30/box)	
	500 mg	8B-S100-HBJ	3 mL (50/box)	
	500 mg	8B-S100-HCH	6 mL (30/box)	
Giga™ Tube				
	500 mg	8B-S100-HDG	12 mL (20/box)	
	1 g	8B-S100-JDG	12 mL (20/box)	
	1 g	8B-S100-JEG	20 mL (20/box)	
	2 g	8B-S100-KEG	20 mL (20/box)	
	5 g	8B-S100-LFF	60 mL (16/box)	

**Tab-less and Teflon® tubes available. Contact Phenomenex for details.



roQ™ Extraction Kits

Wide selection of QuEChERS kits to suit your extraction and clean up needs.

Extraction kits contain fifty easy-pour salt packets and fifty 50 mL stand-alone centrifuge tubes

Description	Unit	Part No.	Price
AOAC 2007.01 Method Extraction Kits			
6.0 g MgSO ₄ , 1.5 g NaOAc	50/pk	KSO-8911*	
EN 15662 Method Extraction Kits			
4.0 g MgSO ₄ , 1.0 g NaCl, 1.0 g SCTD, 0.5 g SCDS	50/pk	KSO-8909*	
Original Non-buffered Method Extraction Kits			
4.0 g MgSO ₄ , 1.0 g NaCl	50/pk	KSO-8910	
6.0 g MgSO ₄ , 1.5 g NaCl	50/pk	KSO-8912	

*AOAC and EN Extraction kits also available in traditional non-collared 50 mL centrifuge tubes, Part No.: KSO-8911-NC and KSO-8909-NC

roQ dSPE Kits

dSPE kits contain pre-weighed sorbents/salts inside 2 mL or 15 mL centrifuge tubes

Description	Unit	Part No.	Price
2 mL dSPE Kits			
150 mg MgSO ₄ , 25 mg PSA, 25 mg C18E	100/pk	KSO-8913	
150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB	100/pk	KSO-8914	
150 mg MgSO ₄ , 25 mg PSA, 7.5 mg GCB	100/pk	KSO-8915	
150 mg MgSO ₄ , 25 mg PSA	100/pk	KSO-8916	
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18E, 50 mg GCB	100/pk	KSO-8917	
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18E	100/pk	KSO-8918	
150 mg MgSO ₄ , 50 mg PSA, 50 mg GCB	100/pk	KSO-8919	
150 mg MgSO ₄ , 50 mg PSA	100/pk	KSO-8920	
15 mL dSPE Kits			
900 mg MgSO ₄ , 150 mg PSA, 150 mg C18E	50/pk	KSO-8921	
900 mg MgSO ₄ , 150 mg PSA, 15 mg GCB	50/pk	KSO-8922	
900 mg MgSO ₄ , 150 mg PSA, 45 mg GCB	50/pk	KSO-8923	
900 mg MgSO ₄ , 150 mg PSA	50/pk	KSO-8924	
1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18E, 400 mg GCB	50/pk	KSO-8925	
1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18E	50/pk	KSO-8926	
1200 mg MgSO ₄ , 400 mg PSA, 400 mg GCB	50/pk	KSO-8927	
1200 mg MgSO ₄ , 400 mg PSA	50/pk	KSO-8928	

MYCOTOXINS

Analysis

Using QUENCHERS, SPE & LC/MS/MS

• Selective

• Rapid

• Repeatable

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions

Trademarks

Kinetex is a registered trademark, Strata-X, Giga, and roQ are trademarks of Phenomenex. API 5500, API 5000, API 4000 are trademarks of AB Sciex Pte. Ltd AB SCIEX™ is being used under license. Barilla is a registered trademark of Barilla G. s E. Fratelli-Societa per Azioni Joint Stock Company Italy.

Disclaimer

Phenomenex is in no way affiliated with Barilla. Comparative separation may not be representative of all applications.

Strata-X is patented by Phenomenex. U. S. Patent No. 7,119,145

© 2015 Phenomenex, Inc. All rights reserved.

BR28511214_L

Explore Your Options for Food Testing

HPLC/UHPLC

- HPLC/UHPLC/PREP
Core-Shell and Fully Porous Products
- Extensive Range of Stationary Phases
- ISO Certified Quality Management

www.phenomenex.com/HPLC



GC and GC/MS

- Wide Selection of Phases
- High Temperature Limits
- Low Bleed Performance
- Enhanced Aqueous Stability

www.phenomenex.com/GC



Sample Prep

- Filtration, SPE, and Vials
- Solutions for All Matrices
- Extensive Method Development Support
- Easily Accessible Web Tools

www.phenomenex.com/SamplePrep



guarantee

If Phenomenex products in this guide do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, send in comparative data within 45 days and keep the product for FREE.

Access the Latest Food Testing Resources



www.phenomenex.com/Food

- Applications
- Technical Notes
- Online Method Development Tools
- On-Demand Webinars
- News and Events